

# Rigel UNI-SiM

# The world's smallest integrated NiBP, SpO2 and patient simulator.

The UNI-SiM is a handheld and battery-operated vital signs simulator capable of undertaking six synchronized vital signs parameters. This enables medical device engineers to quickly, easily and accurately perform NiBP, SpO2, ECG, temperature, IBP and respiration functionality tests simultaneously, using a single portable instrument.

With a fast boot up and single button simulation to repeat the last simulation value in seconds, the UNI-SiM reduces the time taken to test the correct performance of a wide range of medical devices and equipment.

It is easy-to-use and incorporates the full functionality of conventional NiBP and SpO2 simulators with a comprehensive patient simulator.

Compatibility with the Rigel PULS-R universal SpO2 simulation finger creates a truly versatile and valuable tool for every biomed engineer in need of a cost-effective and lightweight solution for testing the performance of vital signs monitors.



# **Key Features**

- Compact and cost-effective
- 6-in-1 vital signs simulation
- Fast start up and single button simulation
- Accurate and real-life simulations
- User definable NiBP simulations
- On-board automation and data storage
- User programmable patient conditions
- Universal SpO2 simulation with PULS-R
- Easy and accurate probe placement with PULS-R

#### **Simulation Functions**

NiBP (systolic and diastolic)

- ECG
- Respiration
- SpO2
- IBP
- Temperature

## **End User Types**

BMETs requiring a cost-effective solution for carrying out performance checks and full PMs on vital signs monitors.

Download your **FREE** introduction to measuring and simulating Vital Signs **www.rigelmedical.com/guides** 

### www.rigelmedical.com/UNI-SiM

**Tel:** +44 (0) 191 587 8730 **Email:** info@rigeImedical.com





## Compact and cost-effective

A highly cost-effective and compact solution for testing 6 of the most common vital signs using a single battery-powered simulator.





## 6-in-1 vital signs simulation

The UNI-SiM incorporates full NiBP simulation including dynamic and static pressure simulation and leak and over-pressure testing with a comprehensive SpO2 and ECG simulator, all in a single hand-held enclosure.

## Fast start up and single button simulation

Automatic power-up of the most recent settings provides simulation of all 6 vital signs with the press of a single button, saving valuable time when setting up the simulator.



www.rigelmedical.com/UNI-SiM

**Tel:** +44 (0) 191 587 8730 **Email:** info@rigeImedical.com







## Accurate and real-life simulations

Fully synchronized simulation signals provide the closest and most accurate representation of a real patient.

## User definable NiBP simulations

User configurable and physiologically correct systolic and diastolic pressures provide a truly universal and accurate NiBP simulator.





## On-board automation and data storage

Simple record management for up to 5,000 asset records. The simulated values from the patient monitor can be entered directly into the UNI-SiM using the built-in keyboard, improving traceability and reducing the need for manual data recording.



www.rigelmedical.com/UNI-SiM

**Tel:** +44 (0) 191 587 8730 **Email:** info@rigeImedical.com



## User programmable patient conditions

Patient specific physiological conditions can be created and stored in the UNI-SiM to provide a highly customisable simulator. A truly versatile tool capable of meeting even the most demanding test protocols.





## Universal SpO2 simulation with PULS-R

Reduce the need for separate accessories with the universal PULS-R SpO2 simulation finger.

This compact SpO2 simulation enables accurate SpO2 simulations in 1% resolution from as low as 30%\* using the pre-programmed manufacturer specific R-curves.\*subject to monitor capability

## **Easy and accurate probe placement with PULS-R** Unique probe placement LED's ensure accurate and correct simulation for each type of SpO2 probe.

The Rigel PULS-R has status LEDs which light up to indicate whether a probe connection has been achieved.





www.rigelmedical.com/UNI-SiM

**Tel:** +44 (0) 191 587 8730 **Email:** info@rigeImedical.com



# **Technical Specifications - UNI-SiM**

#### Non-Invasive Blood Pressure Simulation

Waveform Pulse Volume Heart Rate Integrated Pump Leak Test Chronometer Digital Manometer Pressure Accuracy Pressure Units Oscillometric High, Medium, Low, Paediatric 20 - 300BPM 0 to 350mmHg user configurable User configurable between 0-350mmHg Configurable up to 999 secs 0 - 410mmHg +/- 0.5% FS mmHg, inHg, kg/cm2, cmH2O, mBar, PSI, in H2O and kPa

## Oxygen Saturation Simulation (PULS-R)

Range Repeatability

Heart Rate

Compatibility

Accuracy

ion (PULS-R) 30 to 100% ± 5%\*\* of reading between 30-59% SPO2 ± 3% of reading between 60-89% SPO2 ± 1% of reading between 90-100% SPO2

5 lead simulation including high level output

on Normal Sinus Rhythm (NSR), ST

Elevation, ST Depression, Myocardial

Accuracy of simulation when used with the corresponding R-curves \*Based on using the same probe and monitor setup

\*\*Note that some monitor types might not be able to display low range sats

30-300BPM\*\*\* ± 1BPM Beijing Choice, Criticare, GE Tuffsat, Masimo, Mindray, Nellcor, Nellcor Oximax, Nihon Kohden, Nonin, Novametrix, Philips / HP

\*\*\*Subject to monitor capability

#### ECG Arrhythmia Simulator

Simulation

Heart Rate Accuracy Amplitudes

Accuracy Connection high-level ECG

#### ST Elevation / Depression

Heart Rate Elevation % Elevation Slope

#### Myocardial Infarction Type

Heart rate

**Tall T** Heart Rate T wave Amplitude Infarction, Tall T 20 – 300 BPM ±1 BPM Lead II : 0.5 – 5 mV (in steps of 0.5 mV). Other leads are proportional to Lead II by the following percentages: Lead II : 60 % Lead III : 40 % V1 : 63 % [Reference LA] V2 : 71 % [Reference LA] V3 : 68 % [Reference LA] V4 : 80 % [Reference LA] V5 : 55 % [Reference LA]

±2% 3.5mm jack plug

20 – 300BPM 7%, 13%, 20% Positive, Negative, Flat

V6:49% [Reference LA]

Ischemia, Injury, Infarction, Inferior Infarction 20 – 300 BPM

80 bpm 0 – 1.2mV (steps of 0.1mV)

### Arrhythmia Waveforms (Atrial)

Simulation5 lead simulationAmplitudes0.5 / 1 / 1.5 / 2 / 2.5 / 3 / 3.5 / 4 / 4.5 / 5 mV.Heart rate (where applicable)20 – 300 BPM

#### Atrial

Sinus Arrhythmia (SA), Missing beat, Atrial Flutter (AFLT), Atrial Fibrillation (AFB), Paroxysmal Atrial Tachycardia (PAT), Junctional Premature Contraction

#### Atrial Conduction

First Degree AV Block, Second Degree AV Block - Mobitz I, Second Degree AV Block - Mobitz II, Third Degree AV Block, Right Bundle Branch Block (RBB), Left Bundle Branch Block (LBB), Left Anterior Hemiblock

#### Ventricular

Premature Ventricular Contraction - Intermittent Premature Ventricular Contraction – Continuous, Bigeminy, Trigeminy, Ventricular Flutter (VFLT), Ventricular Fibrillation Fine (VFBF), Ventricular Fibrillation Coarse (VFBC), Monomorphic Ventricular Tachycardia (MVT), Polymorphic Ventricular Tachycardia (PVT) Right Focal (PVC).

Sine, Square, Triangle

0.1 to 0.9Hz (in steps of 0.1Hz)

#### Performance Waveforms Shape

Rate Amplitude Accuracy Shape Rate Amplitude Accuracy

1 to 100Hz (in steps of 1Hz) Lead II : 0.5 – 5 mV (in steps of 0.5 mV). Other leads as above. 2% Pulse 20 ms pulse duration, 4 second delay Lead II : 1 mV. Other leads as above 2 %

# Pacer Waveforms

Available

QRS Pacer Pulse Amplitude Pacer pulse Polarity Pacer pulse width

#### R Wave Detection

Heart Rate R wave width

## **Temperature Simulation**

Simulation Range Accuracy

Default setting

#### **Respiration Simulation**

Rates Base resistances Accuracy Resistance Variations Accuracy Default Settings Apnoea Simulation synchronous atrial, asynchronous atrial, paver only, ventricular pacer, atrial & ventricular pacer

1mV 0.1 – 2 mV positive, negative 0.1 – 2ms

70 bpm 10 – 120 ms (steps of 10ms)

YSI 400 / 700 Static preset at 25, 33, 37and 41°C ±0.1 °C YSI 400 37°C

5, 10, 15, 30, 60, 120, 180 Breaths per Minute 250, 500, 750, 1000 ohms ±5% 0.1, 0.5, 1.0, 1.5 ohms ±10% 15bpm / 250Ω / 0.1Ω 0 - 60 seconds duration 0 - 300 seconds interval.



## www.rigelmedical.com/UNI-SiM

**Tel:** +44 (0) 191 587 8730 **Email:** info@rigeImedical.com



SEAWARD



# **Technical Specifications - UNI-SiM (Continued)**

5µV / V / mmHg

#### Invasive Blood Pressure Simulation

Channels Static Dynamic Accuracy **Excitation Voltage** Impedance Nominal Simulated Sensitivity

2 channels (Channel 2 set at 50% of Channel 1) 0 to 300mmHg 0-300mmHg for Systolic & Diastolic  $\pm 1 \text{mmHg}$ 2 - 16V 350Ω

## **General Specifications - UNI-SiM**

Operation Battery Charger Supply Battery Life

Memory Capacity

Communication

Size (L x W x D)

**Operating Conditions** 

Storage Environment **Environmental Protection** 

Display

Keypad

Weight

Battery cell, in-situ charge 100-240VAC, 50/60Hz 12VDC centre positive 8 hours standby or a maximum of 200 **NiBP** simulations Approx. 5,000 records via Bluetooth Monochrome, 1/4 VGA full graphics Alpha-numeric <1.5kg, <3.5lbs 270 x 110 x 75mm / 10.5 x 4 x 3" 10-30°C, 0-90% RH - NC -15° - +60°C IP 40

#### Service & Warranty

UNI-SiM comes with a free upgraded 24 month warranty (subject to terms and conditions, available at www.rigelmedical.com/register-product)

#### Standard Accessories (supplied with UNI-SiM) ECG snap-on adaptors

Standard Accessories (supplied with PULS-R)

PULS-R universal SpO2 simulation finger

- Carry case
- NiBP tubing kit
- ECG adaptor module

#### **Optional Accessories**

- IBP connect cables
- NiBP accessories
- Temperature connect cables ECG cables and leads

Quick start guide

Power supply

To find out more, visit www.rigelmedical.com/sim-accessories

## **Specifications - PULS-R**

#### Supported Default R Curves

Beijing Choice	Criticare
GE Tuffsat	Masimo
Mindray	Nellcor
Nellcor Oximax	Nihon Kohden
Nonin	Novametrix
Philips / HP	

Heart Rate Setting

30-300BPM (subject to Monitor Compatibility)

Accuracy of Simulation when used with the Corresponding R Curves		
Resolution	Range	Repeatability*
1% steps	30-59%	±5%**
1% steps	60-89%	±3%
1% steps	90-100%	±1%

#### Part Numbers UNI-SiM PULS-R

370A930 399A910

Rev 2.1. 2019

# www.rigelmedical.com/UNI-SiM

Technical specification subject to change without notice.

Tel: +44 (0) 191 587 8730 Email: info@rigelmedical.com



Quick Start Guide