

Rigel Medical, Seaward Group USA,
6304 Benjamin Road, Suite 506, Tampa, FL 33634

Part of  SEAWARD
GROUP USA

Tel: 813-886-2775 **Fax:** 813-886-2754

Email: sales@seaward-groupusa.com **Web:** rigelmedical.com



Innovating Together

Rigel Uni-Therm Accessory List



Rigel Uni-Therm

Accessory List

Foot Switch Adapter Cables				
Part No.	Description	Compatibility	Image	Specification
367A950	Oscilloscope connection lead	All Manufacturers		Oscilloscope SMB to BNC connection lead
367A951	Foot switch adaptor	ERBE ICC350		4 pin male
367A952	Foot switch adaptor	Valleylab™		4 pin male
367A953	Foot switch adaptor	BOWA ARC 200, 300 and 400		7 pin male
367A955	Foot switch adaptor	ForceTriad™		6 pin male

Part No.	Description	Compatibility	Image	Specification
367A956	Foot switch adaptor	ConMed		4 pin male
367A957	Bipolar foot switch adaptor	Valleylab™, Covidien, Megadyne		3 pin male
367A958	Foot switch adaptor	Olympus PSD/UES 40 and 20		5 pin male
367A959	Bipolar foot switch adaptor	Olympus		3 pin male
367A960	Bipolar foot switch adapter	ConMed		3 pin male
367A948	Covidien two-button test pencil	Valleylab™		Test pencil with female banana jack
367A949	Covidien three-button test pencil	ForceTriad™		Test pencil with female banana jack

Neutral Plate Lead Cables				
Part No.	Description	Compatibility	Image	Specification
800A990	Split neutral plate lead with CQM	Aksa, Martin, Eschmann Bovie (Aaron), Integra, Berchtold, Erbe, Aesculap Valleylab		2 pin female with plastic REM control (0.3m length)
800A991	Single neutral plate lead 6.3 Jack	Valleylab™, Ellman, Integra Intern, Berchtold Intern		6.3 Jack (0.3m length)
800A992	Split neutral plate lead with CQM	Valleylab™		2 pin female and REM control (0.3m length)
800A996	Single neutral plate lead w/o CQM	Valleylab™		2 pin female (0.3m length)

Accessory Cables				
Part No.	Description	Compatibility	Image	Specification
800A993	Foot switch cables	All manufacturers		4 x 1.5 meter cables 2 x 4mm plugs 2 x blue and 2 x yellow
800A994	Active cable	All manufacturers		1 meter red cable 2 x 4mm plugs
800A995	Bipolar cable	All manufacturers		1 meter black cable 2 x 4mm plugs
800A997	Accessory cable set	All manufacturers		Consists of: 800A996, 800A995, 800A994, 800A993, 800A992.
286A954	Crocodile clip	All manufacturers		Additional accessory which may be required during performance testing
367A954	Jumper links (per 3)	Rigel Uni-Therm		3 x 4mm jumper links for side panel of Rigel Uni-Therm
PC-05106	NEMA 5-15 to C13 hospital grade power cord	Rigel Uni-Therm		1.8 meter cable, 10amp

Generator Kits			
Part No.	Description	Compatibility	Specification
UTK-CM	ConMed foot switch adapter and cable kit	Sabre 2400, 5000, 7000	Consists of: 367A956, 367A960, 800A997
UTK-EB	ERBE ICC foot switch adapter and cable kit	ERBE ICC*	Consists of: 367A951, 800A997 * In addition to the ERBE ICC kit, an adapter needs to be purchased directly from ERBE (ERBE foot switch adapter from VIO D/S to ICC two-pedal foot switch - Part Number: 20140-004)
UTK-VL1	Covidien / Valleylab foot switch adapter and cable kit	F, FX	Consists of: 367A952, 367A957, 800A997
UTK-VL2	Covidien / Valleylab foot switch adapter and cable kit	ForceTriad™ FT10	Consists of: 367A952, 367A955, 367A957, 800A997
UTK-VL3	Covidien / Valleylab foot switch adapter, cable and test pencil kit	F, FX	Consists of: 367A952, 367A957, 800A997, 367A948
UTK-VL4	Covidien / Valleylab foot switch adapter, cable and test pencil kit	ForceTriad™ FT10	Consists of: 367A952, 367A955, 367A957, 367A948, 367A949, 800A997

Software	
Part No.	Description
383A910	Med-eBase V2 Asset management software
383A914	Med-eBase V2 additional user license
383A915	Med-eBase V2 additional product license (Rigel Uni-Therm)



Data Communication			
Part No.	Description	Image	Specification
396A976	USB download cable		USB cable for data communication with PC and performance analyzers

Printers and Scanners			
Part No.	Description	Image	Specification
339A923	Bluetooth barcode scanner		Bluetooth scanner for asset detail transfer to Rigel products

Carry Case Solutions			
Part No.	Description	Image	Specification
370A960	Rigel Med-eKit Pro		Protective travel case

At Rigel, we are always developing our range of accessories. If you require an accessory which is not on the list, please contact us.